

Course project

Main target of this project was to analyze and visualize eye-tracking data. To achieve that I wrote a java program:

Compile:

```
javac FixationDetection.java
```

Run:

```
java FixationDetection trains.csv subject_ids+ [Options]
```

Options:

- show=subject_number - shows eye path of subject with fixations.
- dt=? default: --dt=70 - sets dispersion threshold
- min=? default: --min=100 - sets minimum fixation duration
- s=? default: --s=90 - sets minimum saccade amplitude
- tolerance=? default: --tolerance=3 - sets noise tolerance

Examples:

```
java FixationDetection trains.csv s8 s18 s28 s4 s14 s24  
default use for group 8
```

```
java FixationDetection trains.csv s8 s18 s28 s4 s14 s24 --dt=50 --min=100  
sets dispersion threshold=50 and minimum duration=100
```

```
java FixationDetection trains.csv s8 --show=0  
shows visualization of first "s8" sample
```

Outputs:

The program always produces 3 files:

means.csv – statistics csv file with following structure: subject_id MFD_true MFD_SD_true MFD_false MFD_SD_false MSA_true MSA_SD_true MSA_false MSA_SD_false MFD_overall MFD_overall_SD MSA_overall MSA_overall_SD

fixations.csv – list of all fixations with following structure: subject_id,known,duration,x,y

saccades.csv – list of all saccades with following structure: subject_id,known,amplitude

The following files: means.csv, fixations.csv, saccades.csv were produced with following parameters:

```
java FixationDetection trains.csv s8 s18 s28 s4 s14 s24
```

dispersion threshold = 50 (50 pixels of dispersion allowed)

minimum duration=100 (100ms if sample rate is 1000 Hz)

tolerance = 3 (up to 2 samples can be omitted in a row after exceeding threshold)

saccade amplitude = 90 (minimum 90 pixels of saccade amplitude)

Algorithm:

Saccade detection – velocity based algorithm

Fixation detection – dispersion threshold algorithm

```

ArrayList<Fixation> dispersionThresholdAlgorithm(String[] subject, double threshold, int
duration, int tolerance_threshold) {
    boolean reset = true;
    int count = 0;
    double minX = 0;
    double maxX = 0;
    double minY = 0;
    double maxY = 0;
    double oldMinX = 0;
    double oldMaxX = 0;
    double oldMinY = 0;
    double oldMaxY = 0;
    int tolerance = 0;
    ArrayList<Fixation> fixations = new ArrayList<>();
    for (int i = 4; i < subject.length; i = i + 2) {
        count++;
        double x = Double.parseDouble(subject[i]);
        double y = Double.parseDouble(subject[i + 1]);
        if (reset) {
            minX = x;
            maxX = x;
            minY = y;
            maxY = y;
            oldMinX = x;
            oldMaxX = x;
            oldMinY = y;
            oldMaxY = y;
            reset = false;
        }
        if (x < minX) {
            oldMinX = minX;
            minX = x;
        }
        if (x > maxX) {
            oldMaxX = maxX;
            maxX = x;
        }
        if (y < minY) {
            oldMinY = minY;
            minY = y;
        }
        if (y > maxY) {
            oldMaxY = maxY;
            maxY = y;
        }
        double dispersion = Math.max(maxX - minX, maxY - minY); //very simplified
        if (dispersion > threshold) {
            tolerance++;
            minX = oldMinX;
            maxX = oldMaxX;
            minY = oldMinY;
            maxY = oldMaxY;
            if (count >= duration && tolerance >= tolerance_threshold) {
                Point2D center = new Point();
                center.setLocation((maxX + minX) / 2, (maxY + minY) / 2);
                fixations.add(new Fixation(center, count, (maxX - minX > maxY - minY ?
maxX - minX : maxY - minY)));
                count = 0;
            }
        }
    }
}

```

```

        reset = true;
    } else if (tolerance >= tolerance_threshold) {
        i = i - (count * 2) + 2;
        count = 0;
        reset = true;
    }
} else {
    tolerance = 0;
}
}
if (count >= duration) {
    Point2D center = new Point();
    center.setLocation((maxX + minX) / 2, (maxY + minY) / 2);
    fixations.add(new Fixation(center, count, (maxX - minX > maxY - minY ? maxX -
minX : maxY - minY)));
}
return fixations;
}

```



```

ArrayList<Saccade> saccadeDetection(String[] subject, double threshold, double
amplitude_threshold, int tolerance_threshold) {
    int count = 0;
    double realDistance = 0;
    Point2D start = new Point();
    Point2D current = new Point();
    Point2D prev = new Point();
    start.setLocation(Double.parseDouble(subject[2]), Double.parseDouble(subject[3]));
    int tolerance = 0;
    ArrayList<Saccade> saccades = new ArrayList<>();
    for (int i = 4; i < subject.length; i = i + 2) {
        count++;
        current.setLocation(Double.parseDouble(subject[i]), Double.parseDouble(subject[i +
1]));
        prev.setLocation(Double.parseDouble(subject[i - 2]), Double.parseDouble(subject[i -
1]));
        double velocity = current.distance(prev);
        double amplitude = current.distance(start);
        if (velocity < threshold) {
            tolerance++;
            if (amplitude >= amplitude_threshold && tolerance >= tolerance_threshold) {
                saccades.add(new Saccade(start, current, count, realDistance));
                count = 0;
                start.setLocation(current);
                realDistance = 0;
            } else if (tolerance >= tolerance_threshold) {
                count = 0;
                start.setLocation(current);
                realDistance = 0;
            }
        } else {
            tolerance = 0;
            realDistance += velocity;
        }
    }
    if (current.distance(start) >= amplitude_threshold) {
        saccades.add(new Saccade(start, current, count, realDistance));
    }
    return saccades;
}



```

Results:

Saccade amplitude for known=true and known=false



Fixation duration for known=true and known=false:



Some visualization of fixations with eye-path

